

Background

Introduction

Spatial transcriptomics, single-cell RNA sequencing (scRNA-Seq), and single-cell DNA methylation sequencing offer complementary views of cellular biology. By integrating these data types, we aim to achieve an unprecedented understanding of a biological system's complexities. This study demonstrates how combining these approaches yields spatially resolved information about cellular epigenetic states and transcriptomic patterns within an efficient and cost-effective workflow.

Methods

Mouse brain slices of varying thickness (10 μ m and 300 μ m) were collected from three distinct brain regions. The 10 μ m samples were subjected to spatial transcriptomics profiling using Vizgen's 500-gene pan-neuron panel on the MERSCOPE platform. The 300 μ m slices were used for transcriptome-wide RNA profiling with ScaleBio's scRNA-Seq kit and genome-wide DNA methylation sequencing using their scMET kit.

Seurat was employed for processing and analyzing both the ScaleBio scRNA-Seq and Vizgen spatial data. Each dataset underwent normalization with SCTransform and clustering via the default Seurat workflow. Subsequently, the datasets were integrated using Seurat's anchor-based canonical correlation analysis (CCA) approach. In contrast, scMET data was processed through ScaleBio's Seq Suite pipeline and visualized using the Amethyst tool.

Results

All three platforms yielded high-quality data, exhibiting robust cell, gene, and transcript recovery. Vizgen data displayed anticipated spatial patterns of cell-type-specific genes, offering a reliable spatial reference for imputation. Unbiased clusters identified in the ScaleBio scRNA data showed clear spatial resolution upon integration with Vizgen data and projection onto brain slices. The observed and imputed spatial patterns demonstrated high concordance with the Allen Brain Atlas when compared at similar brain regions, further validating the quality of the spatial and single-cell data and their integration.




Moreover, the imputed spatial patterns of differentially expressed genes absent from the Vizgen panel but identified as hallmarks of multiple spatially resolved clusters aligned with Allen Brain Atlas data. This confirmed that accurate spatial patterns can be inferred even for genes exclusively present in the single-cell data.

Leveraging single-cell DNA methylation patterns, we successfully identified diverse cell types that corroborated with both the spatial and transcriptomics datasets. Furthermore, clustering based on CH methylation levels helped distinguish between neuronal and non-neuronal cell types.


Conclusions

This study highlights the power of integrating spatial transcriptomics, single-cell RNA sequencing, and single-cell DNA methylation sequencing. By simultaneously examining epigenetic states and transcriptomic patterns at single-cell resolution within a spatial context, we gain a more comprehensive understanding of the intricate cellular landscape of the mouse brain. This integrated approach holds immense promise for unraveling the complex interplay between gene expression, epigenetic regulation, and spatial organization in various biological systems.


Technology Overview




MERSCOPE Shows Expected Spatial Pattern of Known Genes




Scale scRNA-seq Identifies Distinct & Relevant Cell Types



Integration Shows Distinct Spatial Expression of Single-cell Data



Scale scMET Discerns Single-cell Methylomes & Gene Patterns



Conclusions

- Here we demonstrate the data from adjacent brain tissue slices that were run simultaneously through Vizgen MERSCOPE, Scale Bio scRNA-seq and scDNA methylation profiling yielding high quality spatial data, high throughput transcriptome-wide single-cell RNA data, and genome-wide single-cell DNA methylation data.
- We recovered and identified identical cell types across the three datasets demonstrating good concordance between the three platforms.
- We demonstrate the utility of Amethyst in visualizing complex genome wide single-cell DNA methylation data in identifying and curating cell types in a complex sample type.
- Integration of these three datasets provides comprehensive spatial transcriptomic and epigenomic analyses and offers valuable insights into gene expression organization within complex tissues.

References

Rylaarsdam, Lauren E., et al. "Single-cell DNA methylation analysis tool Amethyst reveals distinct noncanonical methylation patterns in human glial cells." *bioRxiv* (2024): 2024-08.

Experimental Design and Data QC for Integration

